

Exempt Patient Specimens - Packaging Checklist

Packaging

- Primary receptacle
- Secondary container
- Outer container
- Absorbent material

Marking and Labeling

- Address label
- "Exempt Animal Specimens" or "Exempt Human Specimens" marking
- Class 9 Miscellaneous Hazard label (shipping with Dry Ice)
- Air Cargo only label (shipping over 50 mL of specimens by cargo aircraft)

Shipping Documentation

- Air Waybill
- Submission Forms

Exempt Animal Specimens